Welcome to Omni , we're thrilled to have you join the family!

At Conquering Horizons, we're passionate about making a real difference in your life. To ensure Omni is made for you, please can you take some measurements and complete and return the details below.

To understand how Omni can be configured for you, please refer to the handbook on our website (QR Code overleaf) or contact Conquering Horizons to discuss any particular requirements or concerns.

Customer name:

Date of assessment:




Assessor name & role:



Diagnosis / primary condition:

Weight (kg):

Height (cm):

Drive preference: Left ☐ Right ☐

Measurement parameters		Description	Measurement (mm)
Width		Whilst seated on a suitable, firm flat surface, measure the width across the hips.	.....
Depth		Buttock to popliteal fossa (rear of knee). Sit back against a firm surface and measure the distance from the face of the back surface to the popliteal fossa (rear of knee).	.....
Height		Whilst sitting on a firm flat surface, with your feet (with footwear) placed firmly on the floor, measure the distance from the base of the heel to the popliteal fossa (rear of the knee).	.....
Backrest angle		Omni is configured with a nominal setting of 95° to 100° and can be adjusted as described in the handbook.	

Armrest height		Whilst seated on a firm, upright surface, with shoulders relaxed and arms bent, measure the distance between the sitting surface and the underside of the elbow.	.....
Armrest length		Whilst seated against a firm, upright surface, with shoulders relaxed and arms bent measure the distance between the face of the back surface and the lower knuckle of the index finger.	.....
Knee flexion		Omni is configured with a nominal setting of 90°. If an alternative setting is required, please contact Conquering Horizons for further assistance.	
Foot flexion		Omni is configured with a nominal setting of 90°. If an alternative setting is required, please refer to the handbook.	

Footrest	Footrests are available in two lengths. (measured from heel to front edge)	Standard 180mm ☐ Extended 250mm ☐
Lapbelt	Attachment points are avialable for lap belts	Lapbelt fixings ☐
Seat back	The seat back is avaible in 3 different heights: Low Shoulder Head rest 	Low back ☐ Shoulder back ☐ Headrest back ☐
Supplementary	If you have additional requirements, please contact Conquering Horizons for further assistance.	

Conquering Horizons Limited
Unit P2 Creech Business Park, Taunton, Somerset UK TA3 5PX
E: andyr@conqueringhorizons.com
W: www.conqueringhorizons.com

Reference: OM-88-00044 v1
© Conquering Horizons Limited 2026